



BD Hypak™ for Biotech

Glass Pre-fillable Syringe

Cover your bases in biotech





Grounded on decades of BD Hypak™ experience, BD Hypak™ for Biotech has been developed to address the key compatibility challenges for biotech drugs and chronic treatments



A reliable Pre-fillable Syringe
from a trusted partner



More than **75%** of major
injectable biotech drugs
are delivered in BD Pre-fillable Syringe solution¹



The BD Hypak™ for Biotech is a reliable and proven Sterile, Clean and ready to Fill (BD SCF™) syringe barrel container featuring key compatibility factors for biotech drug delivery:



Minimized reject rates on customer process through fewer visual and cosmetic defects³

Improved drug compatibility through specified **low tungsten residues levels**^{2,3}

The BD Hypak™ for Biotech is a reliable and proven Sterile, Clean and ready to Fill (BD SCF™) syringe barrel container featuring key compatibility factors for biotech drug delivery:



Controlled **Length Under Flange** and **gliding specifications** foundational for compatibility with secondary devices^{3,4}



BD Physioject™
Disposable Autoinjector



BD UltraSafe Passive™
Needle guard

BD Hypak™ for Biotech offering:

BD SCF™*	1 mL Long
Flange	Cut or Small Round
Needle Cap type	Rigid Needle Shield
Needle diameter	27G – 29G
Needle Wall thickness	Regular, Thin or Special thin
Needle length	1/2inch
Needle type	5 bevels

References and notes

* SCF: Sterile, Clean and ready to Fill (BD SCF™).

1. Internal analysis based on revenues projections in the scope of 15 therapeutic classes in the chronic conditions, BD internal sales (Vendavo 2017) and BD internal market model (InThought 2017).
2. Validation Summary Low Tungsten BD Hypak™ for Biotech Prefillable Syringe Barrels [Internal Study] Pont-de-Claix, FR: Becton, Dickinson and Company; 2012.
3. BD Hypak™ for Biotech customer specifications.
4. BD Hypak™ for Biotech 1mL Long Prefillage Syringe [Internal Study] Pont-de-Claix, FR: Becton, Dickinson and Company; 2013.

BD Medical

Pharmaceutical Systems

United States

1 Becton Drive

Franklin Lakes, NJ 07417

Phone: +1 800 225 3310

Europe

11 rue Aristide-Bergès

38800 Le Pont-de-Claix

France

Phone: +33 4 76 68 36 36 - Fax: +33 4 76 68 35 05

Becton Dickinson France S.A.S - Share capital: 64 719 915 €

RCS Grenoble B 056 501 711

bd.com

BD, the BD Logo, BD Hypak, BD Physioject, BD SCF and BD UltraSafe Passive are trademarks of Becton, Dickinson and Company or its affiliates. © 2020 BD. All rights reserved.

