

BD tools for anticoagulant combination product developers

# BD Preventis™ Needle Shielding System

## Product

BD Preventis™ Plunger rod - available in different colors for customization



Compatible with 0.5 mL and 1 mL BD Hypak™ SCF™ Glass Prefillable Syringes<sup>2</sup>

## Experience

- Established on the market since 2003
- More than 1.5 billion units sold<sup>6</sup>

## Availability

- Samples available on demand
- Commercial availability to be evaluated against requirements

## Key benefits

- Enables one-handed<sup>3</sup>, controlled activation that automatically shields the needle
- Nurses who received instructions were able to successfully recognize the activation of the needle shielding system<sup>3</sup>
- Comes with technical documentation and support<sup>5</sup>
- **Intended users<sup>4</sup>** : Healthcare Workers (HCWs) & Caregivers
- **Intended use** when assembled with a Glass Prefillable syringe (PFS)<sup>4</sup>: Subcutaneous administration of anticoagulants in both clinical and non-clinical environment

Compatible with BD Hypak™ SCF™ Glass PFS<sup>2</sup>

Volume	0.5mL, 1mL
Needle length	½"
Flange design	Cut Flange
BD Preventis™ Plunger rod	Different colors option

BD does not make any claim with regards to BD Preventis™ compliance to sharps injury requirements and its ability to meet safety legislations requirements



## References

1. DVTR20200072, [Internal report], Pont-de-Claix, FR, Becton Dickinson and Company;
2. SD20190231 rev3.0, [Internal report], Pont-de-Claix, FR, Becton Dickinson and Company;
3. Out of 576 uses, 93.8% of nurses say they activate BD Preventis™ with one hand, BDPS-01-13, [Internal report], Pont-de-Claix, FR, Becton Dickinson and Company, July 2001;
4. RMSR 0069\_2011 Rev 02 [Internal report]. Pont de Claix, FR, Becton, Dickinson and Company;
5. BD provides the operational sequence of use for the product to customer;
6. BD analysis [Internal report]. Pont de Claix, FR, Becton, Dickinson and Company; 2021

### **BD Medical**

Pharmaceutical Systems

#### **United States**

1 Becton Drive

Franklin Lakes, NJ 07417

+1 800 225 3310

#### **Europe**

11 rue Aristide-Bergès

38800 Le Pont-de-Claix

France

Phone: +33 4 76 68 36 36 - Fax: +33 4 76 68 35 05

Becton Dickinson France S.A.S - Share capital: 64 719 915 €

RCS Grenoble B 056 501 711

**bd.com**

